














### ►Easy to use



## subscription



## deep learning



## Seven Features of the Image Analysis Web App IAS

Web Applications: Any OS, Remote Work Made Easy

Support for various microscope formats and metadata: information can be reproduced.

Information can be read and positional and dimensional

Batch analysis of multiple specimens and 3D analysis can be performed, contributing to improved work efficiency.

Full-screen display is available: facilitates multi-person discussions

analysis

Solves problems without time lag when you have trouble with operation, etc.


Support for subscription sales format: Helps reduce initial investment

In addition to the above items, we have various other advantages. For more information, please visit our website or contact us.

## **Practical Application Examples**

The IAS can be used effectively in the following areas. If you would like to know more examples of its use, please contact us.

Please feel free to contact us.









The above items are general-purpose fields. It can be effectively used in many other fields.

Cloud-based web application for remote discussion. The direct observation of the data makes it easy to make decisions.

## **Image Analysis Web App IAS Major Points**

# ► Point1 Automatic analysis with deep learning and machine learning






The image analysis web app IAS is equipped with deep learning and machine learning functions, and is useful for analyzing and measuring various samples in the life science, medical, and industrial fields. Samples that were difficult to recognize in the past can now be easily extracted by making full use of AI functions.

# ► Point2 Multiple format support

With support for a wide range of formats, it is possible to reproduce the conditions of particularly complex microscopic imaging. Deep learning is also supported for intuitive image analysis. A high-security server and multiple authentication allow for secure data management.

3D Image





video



## ► Point3 Display of analysis results



NSEdot plothistogram, and heatmap displays are supported, allowing you to visualize the results of high-throughput analysis of multiple samples.



dot plot





histograph

heat map

### [Option]

Remote control options such as automatic microscopes and cameras are also available. With its function, it can be effectively used for equipment operation of collaborators far away and for closed facilities such as clean rooms.



#### ■ List of typical functions

Software Type

Supported container information

display format

Main filter functions Main measurement

Main statistical analysis display functions

Supported File Formats

Supported Formats

web application

Slides, Deschutes, Well Plates, Wafers

Channel 2D 3D Time-lanse

Stitching, shading correction, white and black balance, binarization, focus composition, homogenization, separation, emphasis, super resolution, deconvolution

Number of pieces, luminance, area, color tone, distance, class

NSE, Dot Plot, Histogram, Heat Map

aim, .arf, .psd, .tif, .tiff, .al3d, .gel, .am, .amiramesh, .grey, .grey, .hx, .labels, .img, .hdr, .sif, .png, .afi, .svs, .avi, .exp, .xm, .h5, .1sc, .pic, .raw, .scn, .ims, .img, .cr2, .crw, .jpg, .thm, .wav, .ch5, .vsi, .ets, .pnl, .htd, .log, .c01, .dib, .cxd, .dic, .dcm, .dicom, .jp2, .j2ki, .j2kr, .raw, . .ima, .cr2, .crw, .wav, .dv, .r3d, .r3d, r3d, .dv.log, .r3d.log, .v, .eps, .epsi, .ps, .flex, .mea, .res, .inf, .dm2, .dm3, .dm4, .gif, .naf, .his, .ndpi, . .vms, .txt, .i2i, .hed, .mod, .inr, .ipl, .ipm, .ff, .ics, .ids, .ipw, .xdce, .xlog, .frm, .hdr, .dat, .par, .jpeq, .jpe, .jp2, .j2k, .jpf, .jpk, .jpx, .klb, .xv, .bip, .sxm, .fli, .lim, .msr, .lei, .lir, .scn, .mnc, .mrc, .st, .ali, .map, .rec, .mrcs, .stk, .nd, .scan, .vff, .mrv, .stp, .mng, .nii, .hdr, .nii. gz,. .nrrd..nhdr..nd2..nef..obf..msr..ome.tiff..ome.tif2..ome.tif8..ome.tif8..ome.btf..companion.ome..ome.ome.xml..apl..tnb.. mtb.. .oib, .oif, .pty, .lut, .oir, .dat, .sld, .spl, .liff, .top, pcoraw, .rec, .pcx, .pict, .df3, .im3, .ano, .cfg, .csv, .htm, .rec, .tim, .zpo, .qptiff, .bin, .pbm, .pgm, .ppm, .cfg, .env, .spe, .afm, .mov, .rcpnl, .sm2, .sm3, .spc, .set, sdt, .spi, .xqd, .xqf, .fake, .tif2, .tf8, .btf, .csv, .vws, .pst, .inf, .tfr, .ffr, .zfr, .zfp, .2fl, .tga, .pr3, .dti, .fdf, .bif, .xys, .html, .mvd2, .aisf, .aix, .atsf, .acff,.wat, .bmp, .wiz, .wpi, .czi, .lms, .lsm, .mdb, .zvi, .zip

3i SlideBook, Andor Bio-imaging Division (ABD) TIFF, AIM, Aicona3D, Amershan Bioscience Gel, Amira Mesh, Amnis FlowSight, Analyze 7.5, Andor SIF Animated PNG. Aperio AFI. Aperio SVS TIFF. Applied Precision CellWorX. AVI (Audio Video Interlearve), Axon Raw Format. BD Pathway, Becker & Hickl SPC FIFO. Becker & Hickl SPCImage, Big Date Viewer, Bio-Rad Gel, Bio-Rad PIC, Bio-Rad SCN, Bitplane Imaris, Bruker MRI, Burleigh, Canon DNG, CellH5, Cellomics, cellSens VSI Cellomics, cellSens VSI, CellVoyager, CV7000, DeltaVision, DICOM, ECAT7, EPS (Encapsulated PostScript), Evotec/PerkinElmer OperaFlex, EVOS, FEI FEI TIFF, FITS (Flexible Image Transport System), Gatan Digtal Micrograph, Gatan Digtal Micrograph2, GE MicroCT, GIF (Graphics Interchange Format) Hamamatsu Aguacosmos NAF, Hamamatsu HIS, Hamamatsu ndpi, Hamamatsu VMS, Hitachi S-4800, I2I, ICS (Image Cytometry Standard), Imacon, ImagePro Sequence, ImagePro Workspace, IMAGI, IMOD, Improvision Openlab LIFF, Improvision Openlab Raw, Improvision TIFF, Imspector OBF, InCell 1000/2000 InCell3000, INR, Inveon, Inopath MIBI, IPLab, IPLab, IVision, JEOL, JPEG, JPEG 2000, JPK, JPX, Keller Lab BlockKhoros VIFF (Visualization Image File Format) Bitmap, Kodos Format) BitmapKodak BIPLambert Instruments FILMLaVision ImspectorLeica LCS LEILeica LAS AF LIF (Leica Image File Format), Leica SCN, LEO, Li- Cor L2D, LIM (Laboratory Imaging/Nikon), MetaMorph 7.5 TIFF, MetaMorph Stack (STK), MetaXpress, MIAS (Maia Scientific), Micro-Manager, Mikroscan TIFF, MINIC MRI, Minolta MRW, MNG (Multiple-image Network Graphics), Molecular Imaging, MRC (Medical Research Council), NEF (Nikon Elactronic Format), NIfTI, Nikon Elements TIFF, Nikon EZ-C1 TIFF, Nikon NIS-Elements ND2, NRRD (Nearly Raw Raster Date), Olympus CellR/APL, Olympus FluoView FV1000 Olympus FluoView TIFF, Olympus OIR, Olympus ScanR, Olympus SIS TIFF, OME- TIFF, OME- XML, OMERO Pyramid, Oxford Instruments, PCORAW, PCX (PC Paintbrush ), Perkin Elmer Densitometer, Perkin Elmer Colombus, Perkin Elmer Nuance, Perkin Elmer Operetta, Perkin Elmer UltraVIEW, Perkin Elmer Vectra QPTIFF, Portable Any Map, Adobe Photoshop PSD Portable Any Map, Adobe Photoshop PSD, Photoshop TIFF, PicoQuant Bin, PICT (Macintosh Picture), PNG (Portable Network Graphics), Prairie Technologies TIFF, Princeton Instruments SPE, Quesant, Quick Time Movie, RHK, SBIG, Seiko, SimplePCI &HCImage, SimplePCI &HCImage TIFF, SM Camera, SPIDER, Targa, Text, TIFF (Tagged Image File Format), TillPhotonics TillVision, Topometrix, Trestle, UBM, Unisoku, Varian FDF, Veeco AFM Ventana BIF, VG SAM, VisiTech XYS, Volocity, Volocity Library Clipping, WA-TOP, Windows Bitmap, Woolz, Zeiss Axio CSM, Zeiss AxioVision TIFF, Zeiss AxioVision ZVI (Zeiss Vision Image), Zeiss CZI, Zeiss LSM (Laser Scanning Microscope) 510/710

learning function

Deep learning, machine learning

#### ■ Product name and price

Product Name	Product Number	Unit	Price
Image Analysis Web Application IAS Single Month Contract 500GB	IAS001M1EN	1 license	US\$360
Image Analysis Web Application IAS 3-month contract 500GB	IAS001M3EN	1 license	US\$1,000
Image Analysis Web Application IAS annual subscription 500GB	IAS001Y1EN	1 license	US\$3,600
Image Analysis Web Application IAS annual subscription 2TB	IAS001Y2EN	1 license	US\$14,400
Image Analysis Web Application IAS Workstation	IAS001WSEN	separate meeting	separate meeting

May be used for research purposes only. Cannot be used for diagnostic purposes or procedurally. All trademarks are the property of their respective owners.

For Research Use Only. Not for use in diagnostic procedures. ©2020 Life Analytics Inc. All rights reserved.

All trademarks are the property of Life Analytics and its subsidiaries unless otherwise specified

All trademarks are the property of Life Analytics and its subsidiaries unless otherwise specified. The prices listed are the manufacturer's suggested retail prices as of December 2021.

Consumption tax is not included

Please note that prices, product specifications, appearance, and descriptions are subject to change without notice.



Manufactured by LifeAnalytics

For more information, please contact: support@lifeanalytics.org